



Anti-TMEM62 monoclonal antibody (DCABH-201622)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	TMEM62 (transmembrane protein 62) is a protein-coding gene. Diseases associated with TMEM62 include congenital dyserythropoietic anemia type i, and congenital dyserythropoietic anemia.
Immunogen	A synthetic peptide of human TMEM62 is used for rabbit immunization.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein A
Conjugate	Unconjugated
Applications	WB, ELISA
Size	1 mg
Buffer	In 1x PBS, pH 7.4
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	TMEM62 transmembrane protein 62 [Homo sapiens (human)]
Official Symbol	TMEM62
Synonyms	TMEM62; transmembrane protein 62;
Entrez Gene ID	80021
Protein Refseq	NP_079232
UniProt ID	Q0P6H9

