



Anti-TMEM68 monoclonal antibody (DCABH-201625)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	TMEM68 (transmembrane protein 68) is a protein-coding gene. GO annotations related to this gene include transferase activity, transferring acyl groups.
Immunogen	A synthetic peptide of human TMEM68 is used for rabbit immunization.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein A
Conjugate	Unconjugated
Applications	WB, ELISA
Size	1 mg
Buffer	In 1x PBS, pH 7.4
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	TMEM68 transmembrane protein 68 [Homo sapiens (human)]
Official Symbol	TMEM68
Synonyms	TMEM68; transmembrane protein 68;
Entrez Gene ID	137695
Protein Refseq	NP_001273586
UniProt ID	Q96MH6
Chromosome Location	8q12.1

Function

transferase activity, transferring acyl groups;
