



Magic™ Anti-Staphylococcus type E Staphylococcus Enterotoxin E Monoclonal antibody, Clone N230636 (DCABY-4638)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Specific to Staphylococcal enterotoxin type E
Target	Staphylococcus type E Staphylococcus Enterotoxin E
Immunogen	Staphylococcus Enterotoxin Type E antibody was raised in mouse using Staphylococcal enterotoxin E as the immunogen.
Isotype	IgG2a
Source/Host	Mouse
Species Reactivity	Staphylococcus
Clone	N230636
Purification	Staphylococcus Enterotoxin Type E antibody was purified by Protein A or G affinity purification
Conjugate	Unconjugated
Applications	ELISA Pr* Suggested pair for sandwich ELISA (Capture - Detection): DCABY-4638 - DCABY-4639
Procedure	Matched Antibody Pairs
Format	Liquid
Size	500 µg
Buffer	Supplied as a liquid in PBS, pH 7.4, 3mM KCl, 1.5mM KH ₂ PO ₄ , 140mM NaCl, 8.0mM Na ₂ HPO ₄
Preservative	None
Storage	Store at 4 °C for short term or -20 °C for longer storage. Avoid repeated freeze/thaw cycles.

BACKGROUND

Introduction An enterotoxin (not to be confused with endotoxin) is a protein toxin released by a

microorganism in the intestine. Enterotoxins are chromosomally encoded exotoxins that are produced and secreted from several bacterial organisms. They are often heat-stable, and are of low molecular weight and water-soluble. Enterotoxins are frequently cytotoxic and kill cells by altering the apical membrane permeability of the mucosal (epithelial) cells of the intestinal wall. They are mostly pore-forming toxins (mostly chloride pores), secreted by bacteria, that assemble to form pores in cell membranes. This causes the cells to die.

Keywords

Staphylococcus Enterotoxin Type E; Enterotoxin E; Enterotoxin Type E;
