



**User's Manual**

# Chicken IgY Antigen ELISA Kit



DEIAFY105



2 plate, 5 plate





This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

**Creative Diagnostics**

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)**  **Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)**

---

## PRODUCT INFORMATION

### Intended Use

The Chicken IgY Antigen ELISA Kit is a capture ELISA for the quantitative measurement of Chicken IgY in egg yolk or blood serum.

### Storage

Store all kit components at 2-8°C. Crystal formation may occur in the wash buffer concentrate during prolonged storage at 2-8°C. The crystals can be re-dissolved by swirling the bottle in warm tap water.