



User's Manual

Chrysoidine ELISA Kit

REF DEIANJ25

Σ 96T



RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The Chrysoidine ELISA Kit is an immunoassay for the quantitative and sensitive screening of Chrysoidine in juices and juice drinks.

Reagents And Materials Provided

1. Microtiter plate (12 X 8 strips) coated with antibody, in an resealable aluminum pouch
2. Calibrators/Standards
3. Chrysoidine-HRP Conjugate
4. Anti-Chrysoidine Antibody Solution
5. Wash Solution
6. Substrate (Color) Solution (TMB)
7. Stop Solution
8. Sample Diluent

Storage

Store the kit at 2 - 8°C.