



# Anti-Gentamycin monoclonal antibody, clone 3.DD.906 (DMABT-Z60374)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	Mouse Anti-Gentamycin Monoclonal Antibody
<b>Target</b>	Gentamycin
<b>Immunogen</b>	Gentamycin conjugated to KLH.
<b>Isotype</b>	IgG2a
<b>Source/Host</b>	Mouse
<b>Species Reactivity</b>	N/A
<b>Clone</b>	3.DD.906
<b>Purification</b>	Protein G purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	ELISA
<b>Format</b>	Liquid
<b>Size</b>	100 µg
<b>Buffer</b>	PBS, pH 7.4
<b>Preservative</b>	None
<b>Storage</b>	Store at +4°C short term. Aliquot and store at -20°C or -80°C. Store in the dark. Avoid repeated freeze / thaw cycles.

## BACKGROUND

**Introduction** Gentamycin is an aminoglycoside antibiotic used to treat infections caused by many different types of bacteria. Gentamycin is usually administered by intravenous infusion or intramuscular injection.

**Keywords** Gentamicin;