



Anti-Gentamycin monoclonal antibody, clone 85.5 (DMABT-Z60464)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Mouse Anti-Gentamycin Monoclonal Antibody
Target	Gentamycin
Immunogen	Gentamicin conjugated to KLH.
Isotype	IgG2a
Source/Host	Mouse
Species Reactivity	N/A
Clone	85.5
Purification	Protein G purified
Conjugate	Unconjugated
Applications	ELISA, IA
Format	Liquid
Size	500 µg
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	Gentamycin is an aminoglycoside antibiotic used to treat infections caused by many different types of bacteria. Gentamycin is usually administered by intravenous infusion or intramuscular injection.
Keywords	Gentamicin;