



Anti-DIXDC1A polyclonal antibody (DPAB-DC1712)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide corresponding to zebrafish dixdc1.
Source/Host	Rabbit
Species Reactivity	Zebrafish
Conjugate	Unconjugated
Applications	WB, ICC, IF,
Format	Liquid
Size	25 µg
Buffer	In PBS (0.1% proclin, 2.0% Block Ace)
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	dixdc1a DIX domain containing 1a [Danio rerio (zebrafish)]
Official Symbol	DIXDC1A
Synonyms	DIXDC1A; DIX domain containing 1a; ccd1; dixdc1; dixin-A; dixin; coiled-coil-DIX1; coiled-coil DIX 1; coiled-coil-DIX1-A; DIX domain containing 1; coiled-coil protein DIX1; coiled-coil protein DIX1-A; DIX domain-containing protein 1; DIX domain-containing protein 1-A;
Entrez Gene ID	360137
mRNA Refseq	NM_182884
Protein Refseq	NP_878304
UniProt ID	Q804T6
Chromosome Location	chromosome: 15
Function	gamma-tubulin binding; mitogen-activated protein kinase kinase kinase binding; signal transducer activity;