



Anti-Streptomycin polyclonal antibody (DPAB-DC4776)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Streptomycin conjugated with BSA.
Source/Host	Sheep
Species Reactivity	N/A
Conjugate	Unconjugated
Applications	EIA,
Format	Liquid
Size	100 µl
Buffer	In 20 mM phosphate buffer, 150 mM NaCl, pH 7.2 (0.09% sodium azide)
Preservative	0.09% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Streptomycin is an antibiotic (antimycobacterial) drug, the first of a class of drugs called aminoglycosides to be discovered, and it was the first cure for tuberculosis. It is derived from the actinobacterium <i>Streptomyces griseus</i> . Streptomycin is a bacte
Keywords	Streptomycin