



# Anti-Decoquinatate polyclonal antibody (DPAB27721)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	Sheep Anti-Decoquinatate Polyclonal Antibody
<b>Target</b>	Decoquinatate
<b>Immunogen</b>	Decoquinatate-BSA
<b>Source/Host</b>	Sheep
<b>Species Reactivity</b>	N/A
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	IA, ELISA, WB, Dot, IHC
<b>Format</b>	Liquid/Ig Fraction
<b>Concentration</b>	10.7mg/ml
<b>Size</b>	1 ea
<b>Buffer</b>	20mM Phosphate, 150mM Sodium
<b>Preservative</b>	0.09% Sodium Azide
<b>Storage</b>	Store at 4°C for up to 3 months. For longer term storage aliquot into small volumes and store at -20°C.

## BACKGROUND

<b>Introduction</b>	Decoquinatate is a coccidiostat used in veterinary medicine.
<b>Keywords</b>	6-Decoxy-7-ethoxy-4-oxo-1H-quinoline-3-carboxylic acid ethyl ester; Sheep Anti-Decoquinatate Polyclonal Antibody; Anti-Decoquinatate Polyclonal Antibody; Decoquinatate Polyclonal Antibody Sheep Anti-Decoquinatate PAb; Anti-Decoquinatate PAb; Decoquinatate PAb; Sheep