



Anti-Toltrazuril polyclonal antibody (DPABH-00048)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Sensitivity: 10 ng/ml Toltrazuril produces 53% inhibition in a competitive ELISA, employing Toltrazuril polyclonal antibody.
Immunogen	Toltrazuril conjugated to Bovine Thyroglobulin.
Isotype	IgG
Source/Host	Sheep
Species Reactivity	Bovine
Purity	IgG fraction
Conjugate	Unconjugated
Applications	ELISA
Format	Liquid
Size	50 µl
Buffer	pH: 7.20; Constituents: 99% Phosphate Buffer, 0.88% Sodium chloride
Preservative	0.09% Sodium Azide
Storage	Store at 4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	Toltrazuril is a coccidiostat, an antiprotozoal agent that acts upon Coccidia parasites.
Keywords	Toltrazuril