



Anti-SMYD5 (aa 305-401) polyclonal antibody (DPABH-00660)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	SMYD5 (SMYD family member 5) is a protein-coding gene. GO annotations related to this gene include molecular_function and metal ion binding.
Immunogen	Synthetic peptide, corresponding to a region within C terminal amino acids 305-401 of Human SMYD5.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.00; Constituents: 0.75% Glycine, 1.21% Tris, 10% Glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	SMYD5 SMYD family member 5 [Homo sapiens]
Official Symbol	SMYD5
Synonyms	SMYD5; SMYD family member 5; RAI15, retinoic acid induced 15; SET and MYND domain-containing protein 5; NN8 4AG; RRG1; ZMYND23; retinoic acid induced 15; retinoic acid responsive; retinoic acid-induced protein 15; RAI15; NN8-4AG;

Entrez Gene ID	10322
Protein Refseq	NP_006053
UniProt ID	Q6GMV2
Chromosome Location	2p12
Function	metal ion binding; molecular_function;