



Anti-REEP5 (aa 140-169) polyclonal antibody (DPABH-01436)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May promote functional cell surface expression of olfactory receptors.
Immunogen	Synthetic peptide conjugated to KLH, corresponding to a region within the C terminal amino acids 140-169 of Human REEP5 (NP_005660.4).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	400 µl
Buffer	Constituent: 99% PBS
Preservative	0.09% Sodium Azide
Storage	Store at 4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	REEP5 receptor accessory protein 5 [Homo sapiens]
Official Symbol	REEP5
Synonyms	REEP5; receptor accessory protein 5; C5orf18, chromosome 5 open reading frame 18; receptor expression-enhancing protein 5; D5S346; deleted in polyposis 1; DP1; polyposis coli region hypothetical protein DP1; polyposis locus protein 1; TB2; receptor expression enhancing protein 5; YOP1; C5orf18; MGC70440;

Entrez Gene ID	7905
Protein Refseq	NP_005660
UniProt ID	Q00765
Chromosome Location	5q22-q23
Function	molecular_function; protein binding;