



Anti-TBC1D20 (aa 334-363) polyclonal antibody (DPABH-01509)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May act as a GTPase-activating protein for Rab family protein(s).
Immunogen	A synthetic peptide, conjugated to KLH, corresponding to a region within C terminal amino acids 334-363 of Human TBC1D20 (NP_653229.1).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	400 µl
Buffer	Constituent: 99% PBS
Preservative	0.09% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	TBC1D20 TBC1 domain family, member 20 [Homo sapiens]
Official Symbol	TBC1D20
Synonyms	TBC1D20; TBC1 domain family, member 20; C20orf140, chromosome 20 open reading frame 140; TBC1 domain family member 20; dJ852M4.2; C20orf140; FLJ45119;
Entrez Gene ID	128637
Protein Refseq	NP_653229

UniProt ID	Q96BZ9
Chromosome Location	20p13
Function	GTPase activator activity; Rab GTPase activator activity;