



Rabbit anti-Human MUPCDH polyclonal antibody (DPABH-02775)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Acts as a calcium-dependent cell adhesion protein.
Specificity	BLAST analysis of the peptide immunogen showed no homology with other human proteins.
Immunogen	Synthetic peptide, corresponding to 17 amino acids from the internal region of Human MUPCDH. Q9HBB8.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Format	Liquid
Size	50 µg
Buffer	Constituent: 99% PBS
Preservative	Preservative: 0.1% Sodium azide
Storage	Store at -70°C for long-term. Store at 4°C for short-term.

GENE INFORMATION

Gene Name	CDHR5 cadherin-related family member 6 [Homo sapiens]
Official Symbol	CDHR5
Synonyms	CDHR5; cadherin-related family member 5; MUCDHL; MUPCDH; MU-PCDH; mu-protocadherin; mucin and cadherin-like protein;
Entrez Gene ID	53841

mRNA Refseq	NM_001171968.1
Protein Refseq	NP_001165439.1
MIM	606839
UniProt ID	B4DV98
Function	beta-catenin binding; calcium ion binding;