



Anti-CCDC151 (aa 46-317) polyclonal antibody (DPABH-03504)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	CCDC151 (coiled-coil domain containing 151) is a protein-coding gene.
Immunogen	Recombinant fragment corresponding to a region within amino acids 46-317 of CCDC151 (A5D8V7).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	100 µl
Buffer	pH: 7.00; Constituents: 78% PBS, 20% Glycerol, 1% BSA
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	CCDC151 coiled-coil domain containing 152 [Homo sapiens]
Official Symbol	CCDC151
Synonyms	CCDC151; coiled-coil domain containing 151; coiled-coil domain-containing protein 151;
Entrez Gene ID	115948
Protein Refseq	NP_659482.3
UniProt ID	A5D8V7

Function

protein binding;
