



# Anti-SAP30BP (aa 25-129) polyclonal antibody (DPABH-04075)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Antigen Description</b>	Induces cell death. May act as a transcriptional corepressor of a gene related to cell survival. May be involved in the regulation of beta-2-microglobulin genes.
<b>Immunogen</b>	Synthetic peptide, corresponding to a region within amino acids 25-129 of Human SAP30BP.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Purification</b>	Immunogen affinity purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	ICC/IF, WB
<b>Format</b>	Liquid
<b>Size</b>	100 µl
<b>Buffer</b>	pH: 7.00; Constituents: 0.75% Glycine, 1.21% Tris, 10% Glycerol
<b>Preservative</b>	None
<b>Storage</b>	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

## GENE INFORMATION

<b>Gene Name</b>	<a href="#">SAP30BP SAP31 binding protein [ Homo sapiens ]</a>
<b>Official Symbol</b>	SAP30BP
<b>Synonyms</b>	SAP30BP; SAP30 binding protein; HTRG; HTRP; HCNGP; SAP30-binding protein; HSV-1 binding; transcriptional regulator protein HCNGP;
<b>Entrez Gene ID</b>	<a href="#">29115</a>
<b>Protein Refseq</b>	<a href="#">NP_037392.1</a>

<b>UniProt ID</b>	<a href="#">Q9UHR5</a>
<b>Pathway</b>	Epigenetic regulation of gene expression; Negative epigenetic regulation of rRNA expression;
<b>Function</b>	protein binding;