



Anti-THAP3 (full length) polyclonal antibody (DPABH-06975)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	THAP3 contains 1 THAP-type zinc finger. There are three known isoforms produced by alternative splicing.
Immunogen	Recombinant full length protein, corresponding to amino acids 1-239 of Human THAP3 (AAH92427.1, UniProt ID: Q8WTV1).
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Human
Purification	Protein A purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	50 µg
Buffer	pH: 7.20; Constituent: 99% PBS
Preservative	None
Storage	Store at 4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	THAP3 THAP domain containing, apoptosis associated protein 3 [Homo sapiens]
Official Symbol	THAP3
Synonyms	THAP3; THAP domain containing, apoptosis associated protein 3; THAP domain-containing protein 3; MGC33488;

Entrez Gene ID	90326
Protein Refseq	NP_001182681
UniProt ID	Q8WTV1
Chromosome Location	1p36.1
Function	DNA binding; metal ion binding;