



Rabbit Anti-Human TRIM14 Polyclonal Antibody (DPABH-08000)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant Protein, antigen sequence: HIDNITQIEDATEKCLKANAESSKTWLGKGFTELRLLLDDEEEALAKKFIDKNTQLTLQVYREQADSC REQLDIMNDLS
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen affinity purified
Conjugate	Unconjugated
Applications	IHC, ICC-IF
Format	Liquid
Size	100 µL
Buffer	40% glycerol and PBS (pH 7.2).
Preservative	0.02% Sodium Azide
Storage	Store at +4°C for short term storage. Long time storage is recommended at -20°C. Gently mix before use. Optimal concentrations and conditions for each application should be determined by the user.

BACKGROUND

Introduction	The protein encoded by this gene is a member of the tripartite motif (TRIM) family. The TRIM motif includes three zinc-binding domains, a RING, a B-box type 1 and a B-box type 2, and a coiled-coil region. The protein localizes to cytoplasmic bodies and its function has not been determined. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Mar 2010]
---------------------	---

Keywords

TRIM14; tripartite motif containing 14; tripartite motif-containing protein 14; tripartite motif-containing 14; tripartite motif protein TRIM14;

GENE INFORMATION

Entrez Gene ID [9830](#)

UniProt ID [Q14142](#)
