



# Anti-PCED1B (aa 239-288) polyclonal antibody (DPABH-12652)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Antigen Description</b>	PCED1B (PC-esterase domain containing 1B) is a protein-coding gene. GO annotations related to this gene include hydrolase activity. An important paralog of this gene is PCED1A.
<b>Immunogen</b>	Synthetic peptide within Human FAM113B aa 239-288 (C terminal). The exact sequence is proprietary. (NP_612380)Sequence: CLSQLLLAHVADAWGVLPVGEWIKKKKPGPRVEGPPQANRNHPAL Database link: Q96HM7
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Purification</b>	Immunogen affinity purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB
<b>Format</b>	Liquid
<b>Size</b>	50 µg
<b>Buffer</b>	Constituents: 98% PBS, 2% Sucrose
<b>Preservative</b>	None
<b>Storage</b>	Shipped at 4°C. Store at 4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.

## GENE INFORMATION

<b>Gene Name</b>	<a href="#">PCED1B PC-esterase domain containing 2B [ Homo sapiens ]</a>
<b>Official Symbol</b>	PCED1B
<b>Synonyms</b>	PCED1B; PC-esterase domain containing 1B; FAM113B; PC-esterase domain-containing protein

1B; family with sequence similarity 113, member B;

---

<b>Entrez Gene ID</b>	<a href="#">91523</a>
<b>Protein Refseq</b>	<a href="#">NP_001268358.1</a>
<b>UniProt ID</b>	<a href="#">A0A024R115</a>
<b>Function</b>	hydrolase activity;

---