



Anti-PTRHD1 (aa 76-140) polyclonal antibody (DPABH-19313)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|----------------------------|---|
| Antigen Description | PTRHD1 (peptidyl-tRNA hydrolase domain containing 1) is a protein-coding gene. GO annotations related to this gene include aminoacyl-tRNA hydrolase activity and protein tyrosine phosphatase activity. |
| Immunogen | Recombinant fragment corresponding to amino acids 76-140 (GRMRKVVLEA PDETTLKELA ETLQQKNIDH MLWLEQPENI ATCIALRPYP KEEVGQYLKK FRLFK) of Human PTRHD1 (Q6GMV3). |
| Isotype | IgG |
| Source/Host | Rabbit |
| Species Reactivity | Human |
| Purification | Immunogen affinity purified |
| Conjugate | Unconjugated |
| Applications | ICC/IF, IHC-P |
| Format | Liquid |
| Size | 100 µl |
| Buffer | pH: 7.20; Constituents: 59% PBS, 40% Glycerol |
| Preservative | 0.02% Sodium Azide |
| Storage | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. |

GENE INFORMATION

| | |
|------------------------|---|
| Gene Name | PTRHD1 peptidyl-tRNA hydrolase domain containing 2 [Homo sapiens] |
| Official Symbol | PTRHD1 |
| Synonyms | PTRHD1; peptidyl-tRNA hydrolase domain containing 1; C2orf79; putative peptidyl-tRNA hydrolase PTRHD1; peptidyl-tRNA hydrolase domain-containing protein 1; |

| | |
|-----------------------|---|
| Entrez Gene ID | 391356 |
| Protein Refseq | NP_001013685.1 |
| UniProt ID | Q6GMV3 |
| Function | aminoacyl-tRNA hydrolase activity; protein tyrosine phosphatase activity; |