



Anti-SAP30BP (aa 90-140) polyclonal antibody (DPABH-20000)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Induces cell death. May act as a transcriptional corepressor of a gene related to cell survival. May be involved in the regulation of beta-2-microglobulin genes.
Immunogen	Synthetic peptide, corresponding to a region within amino acids 90-140 of Human SAP30BP (NP_037392.1)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, IP
Format	Liquid
Size	100 µg
Buffer	Constituents: Tris citrate/phosphate, pH 7-8
Preservative	Preservative: 0.09% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	SAP30BP SAP31 binding protein [Homo sapiens]
Official Symbol	SAP30BP
Synonyms	SAP30BP; SAP30 binding protein; HTRG; HTRP; HCNGP; SAP30-binding protein; HSV-1 binding; transcriptional regulator protein HCNGP;
Entrez Gene ID	29115

mRNA Refseq	NM_013260.6
Protein Refseq	NP_037392.1
MIM	610218
UniProt ID	Q9UHR5
Pathway	Epigenetic regulation of gene expression; Negative epigenetic regulation of rRNA expression;
Function	protein binding;