



Anti-SPATA31E1 (aa 913-1032) polyclonal antibody (DPABH-20121)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	SPATA31E1 (SPATA31 subfamily E, member 1) is a protein-coding gene. An important paralog of this gene is SPATA31A1.
Immunogen	antigen, corresponding to amino acids 913-1032 of Human C9orf79 (Q6ZUB1).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	SPATA31E1 SPATA31 subfamily E, member 2 [Homo sapiens]
Official Symbol	SPATA31E1
Synonyms	SPATA31E1; SPATA31 subfamily E, member 1; C9orf79; FAM75E1; spermatogenesis-associated protein 31E1; XXyac-YM21GA2.5; protein FAM75E1; FAM75-like protein C9orf79; family with sequence similarity 75, member E1;
Entrez Gene ID	286234

Protein Refseq [NP_849150.3](#)

UniProt ID [Q6ZUB1](#)
