



# Rabbit Anti-Human THSD4 Polyclonal antibody (DPABH-20213)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	THSD4 Fusion Protein expressed in E. coli
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Purification</b>	Antigen affinity purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	IHC, ELISA Recommended concentration: IHC: 1:50-1:500
<b>Molecular Weight</b>	55 kDa
<b>Positive Control</b>	Transfected HEK-293 cells
<b>Format</b>	Liquid
<b>Size</b>	150 µl
<b>Buffer</b>	PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
<b>Preservative</b>	0.02% sodium azide
<b>Storage</b>	Store at -20°C. Aliquoting is unnecessary for -20°C storage.

## BACKGROUND

**Keywords** THSD4; thrombospondin, type I, domain containing 4; ADAMTSL6; FVSY9334; PRO34005; ADAMTSL-6; thrombospondin type-1 domain-containing protein 4; ADAMTS-like protein 6; A disintegrin and metalloproteinase with thrombospondin motifs-like protein 6;

# GENE INFORMATION

Entrez Gene ID [79875](#)

---

UniProt ID [Q6ZMP0](#)

---