



Anti-FAM65B (aa 765-863) polyclonal antibody (DPABH-20459)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Isoform 2 play a role in promoting myogenic cell differentiation, cytoskeletal rearrangement and filopodia formation.
Immunogen	antigen sequence: TPFVARSLLK KLSRQIQVME KLAADVSDENI GNISVVVEAI PEFHKKLSLL SFWTKCCSPV GYHSPADRV MKQLEASFAR TVNKEYPGLA DPVFRTLVS, corresponding to amino acids 765-863 of Human FAM65B.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	FAM65B family with sequence similarity 65, member B [Homo sapiens]
Official Symbol	FAM65B
Synonyms	FAM65B; family with sequence similarity 65, member B; C6orf32, chromosome 6 open reading frame 32; protein FAM65B; DIFF48; KIAA0386; PL48; DIFF40; C6orf32;

Entrez Gene ID	9750
Protein Refseq	NP_055537
UniProt ID	B7Z5J9
Chromosome Location	6p22.3-p21.32
Function	binding;