



# Anti-CCDC152 (aa 158-254) polyclonal antibody (DPABH-20619)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Antigen Description</b>	CCDC152 (coiled-coil domain containing 152) is a protein-coding gene.
<b>Immunogen</b>	antigen: KEMEISELNA KLRSQEKEKQ NEIIKLQLEF DAKLARVQTK SKSYQDSTVL PQSIYRRKLQ HFQEEKNKEI AILRNTIRDL EQRLSVGKDS HLKRRRF, corresponding to amino acids 158-254 of Human CCDC152 (Q4G0S7).
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Purification</b>	Immunogen affinity purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	ICC/IF, IHC-P
<b>Format</b>	Liquid
<b>Size</b>	100 µl
<b>Buffer</b>	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
<b>Preservative</b>	0.02% Sodium Azide
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

## GENE INFORMATION

<b>Gene Name</b>	<a href="#">CCDC152 coiled-coil domain containing 153 [ Homo sapiens ]</a>
<b>Official Symbol</b>	CCDC152
<b>Synonyms</b>	CCDC152; coiled-coil domain containing 152; CH5400; coiled-coil domain-containing protein 152;
<b>Entrez Gene ID</b>	<a href="#">100129792</a>

Protein Refseq [NP\\_001128320.1](#)

---

UniProt ID [A0A024R043](#)

---