



Anti-PP2D1 (aa 487-577) polyclonal antibody (DPABH-20674)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	antigen sequence: NKSPSKGPLL FSTSEPNLTK SQSNIHVLQ YKSVSEVRVS TTNSKENLSD SNYSKYCIYN PENVETFP AE TTHRKCSEK VTDRPTSVND V, corresponding to amino acids 487-577 of Human PP2D1 (UniProt ID: A8MPX8).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P, WB, ICC/IF
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	PP2D1 protein phosphatase 2C-like domain containing 2 [Homo sapiens]
Official Symbol	PP2D1
Synonyms	PP2D1; protein phosphatase 2C-like domain containing 1; C3orf48; protein phosphatase 2C-like domain-containing protein 1;
Entrez Gene ID	151649
Protein Refseq	NP_001239586.1

UniProt ID

[A8MPX8](#)

Function

catalytic activity;
