



Anti-TBC1D12 (aa 264-371) polyclonal antibody (DPABH-20756)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May act as a GTPase-activating protein for Rab family protein(s).
Immunogen	antigen sequence, corresponding to amino acids 264-371 of Human TBC1D12.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P, ICC/IF
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	TBC1D12 TBC1 domain family, member 13 [Homo sapiens]
Official Symbol	TBC1D12
Synonyms	TBC1D12; TBC1 domain family, member 12; TBC1 domain family member 12;
Entrez Gene ID	23232
Protein Refseq	NP_056003.1
UniProt ID	O60347

Function

Rab GTPase activator activity;
