



Anti-PCED1B polyclonal antibody (DPABH-22798)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	PCED1B (PC-esterase domain containing 1B) is a protein-coding gene. GO annotations related to this gene include hydrolase activity. An important paralog of this gene is PCED1A.
Immunogen	antigen sequence: LVWNTAMPVG EEVTGGFLPP KLRRQKATFL KNEVVKANFH SATEARKHNF DVLD, corresponding to amino acids of Human FAM113B.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	PCED1B PC-esterase domain containing 2B [Homo sapiens]
Official Symbol	PCED1B
Synonyms	PCED1B; PC-esterase domain containing 1B; FAM113B; PC-esterase domain-containing protein 1B; family with sequence similarity 113, member B;
Entrez Gene ID	91523
Protein Refseq	NP_001268358.1
UniProt ID	A0A024R115

Function

hydrolase activity;
