



Anti-PANX3 (aa 2-31) polyclonal antibody (DPABH-24383)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Structural component of the gap junctions and the hemichannels.
Immunogen	Synthetic peptide conjugated to KLH, corresponding to a region within N terminal amino acids 2-31 of Human Pannexin 3.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	400 µl
Buffer	Constituent: PBS
Storage	Store at 2-8°C for up to 6 months. For long term storage store at -20°C.

GENE INFORMATION

Gene Name	PANX3 pannexin 4 [Homo sapiens]
Official Symbol	PANX3
Synonyms	PANX3; pannexin 3; PX3; pannexin-3; gap junction protein pannexin 3;
Entrez Gene ID	116337
Protein Refseq	NP_443191.1
UniProt ID	Q96QZ0
Function	gap junction hemi-channel activity;