



Anti-Nandrolone polyclonal antibody (DPABY-783)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|---------------------------|--|
| Immunogen | 19 Nortestosterone 17 antibody was raised in sheep using 19 Nortestosterone(17)-BSA as the immunogen |
| Source/Host | Sheep |
| Species Reactivity | N/A |
| Conjugate | Unconjugated |
| Applications | ELISA Pr* |
| Format | Liquid |
| Size | 100 µl |
| Buffer | Ig fraction supplied in liquid form in 20 mM phosphate buffer, 150 mM sodium chloride and 0.09% sodium azide, pH 7.2 |
| Preservative | 0.09% Sodium Azide |
| Storage | Store at 2-8 °C for short term storage or -20 °C for up to 2 years |

BACKGROUND

Introduction Nortestosterone, also known as nandrolone, belongs to the group of anabolic androgenic steroids. In some cases 19-nortestosterone is a naturally occurring compound, especially in male pigs, and pregnant cows. The EU banned the use of growth promoters in meat production since 1986.

Keywords 19 Nortestosterone 17; 17b 19-nortestosterone; Nandrolone; Nortestosterone;