



Anti-Neomycin polyclonal antibody (DPABY-922)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|---------------------------|--|
| Immunogen | Neomycin antibody was raised in sheep using Neomycin-BTG as the immunogen |
| Source/Host | Sheep |
| Species Reactivity | N/A |
| Conjugate | Unconjugated |
| Applications | ELISA Pr* |
| Format | Liquid |
| Size | 100 µl |
| Buffer | Ig fraction supplied in liquid form in 20 mM phosphate buffer, 150 mM sodium chloride and 0.09% sodium azide, pH 7.2 |
| Preservative | 0.09% Sodium Azide |
| Storage | Store at 2-8 °C for short term storage or -20 °C for up to 2 years |

BACKGROUND

Introduction Neomycin is overwhelmingly used as a topical preparation, such as Neosporin. It can also be given orally, where it is usually combined with other antibiotics.

Keywords Neomycin; Neomycin Sulfate;