



Anti-Spectinomycin polyclonal antibody (DPABY-964)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Spectinomycin antibody was raised in sheep using Spectinomycin-BSA as the immunogen
Source/Host	Sheep
Species Reactivity	N/A
Conjugate	Unconjugated
Applications	ELISA Pr*
Format	Liquid
Size	100 µl
Buffer	Ig fraction supplied in liquid form in 20 mM phosphate buffer, 150 mM sodium chloride and 0.09% sodium azide, pH 7.2
Preservative	0.09% Sodium Azide
Storage	Store at 2-8 °C for short term storage or -20 °C for up to 2 years

BACKGROUND

Introduction	Aminoglycosides are a family of bacterial antibiotics that are used in the treatment of specific bacterial infections. They display a concentration dependent killing action and are active against a wide range of aerobic Gram-negative bacilli.
Keywords	Spectinomycin; Spectinomycin hydrochloride; Trobicin;