



Magic™ Anti-C21ORF78.6 (full length) polyclonal antibody (DPATB-H81279)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Chicken Anti-C21orf78.6 Polyclonal Antibody
Target	C21orf78.6
Immunogen	Recombinant full length C21orf78.6(Human)
Isotype	IgY
Source/Host	Chicken
Species Reactivity	Human
Purification	Protein L purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	100 µg
Buffer	PBS
Preservative	0.02% Sodium Azide
Storage	Store at +4°C short term. Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	C21orf78.6 is a novel 40 kDa protein encoded on human chromosome 21.
Keywords	C21orf78.6