



Anti-Diclazuril polyclonal antibody (DPATB-H81377)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Sheep Anti-Diclazuril Polyclonal Antibody
Target	Diclazuril
Immunogen	Diclazuril (4-N)-BTG
Isotype	IgG
Source/Host	Sheep
Species Reactivity	N/A
Purification	IgG fraction
Conjugate	Unconjugated
Applications	EIA
Format	Liquid
Size	50 µl
Buffer	20mM PBS, 150mM Sodium chloride, pH 7.2
Preservative	0.09% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	Diclazuril is an anticoccidial drug used in major poultry species, such as broiler chickens, replacement pullets and turkeys, as well as in rabbits, lambs and minor edible bird species. It has not previously been reviewed by the Committee. The chemical na
Keywords	Diclazuril; Nuoqui; Vecoxan; Diclazuril premix 0.5%,1%; CLINACOX; DICLAZURIL