



Anti-E. coli O125 Polyclonal antibody (DPATB-H81424)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-E. coli O125 Polyclonal Antibody
Target	E. coli O125
Immunogen	Tissue / cell preparation Escherichia coli O serotype 125 (heat inactivated).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	E. coli
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	Agg
Positive Control	Antigenic suspension.
Format	Liquid
Size	2 ml
Buffer	Whole serum
Preservative	None
Storage	Store at +4°C. Do not freeze.

BACKGROUND

Introduction	Escherichia coli is a gram negative bacillus that belongs to a larger group of Enterobacteriaceae bacteria that inhabit the gastrointestinal tract. Although usually a harmless resident of the gut, some strains have the potential to cause serious problems, es
Keywords	E. coli O125; Bacteria; Proteobacteria; Gammaproteobacteria; Enterobacteriales; Enterobacteriaceae; Escherichia; E. coli; Escherichia coli; B