



Anti-Tobramycin polyclonal antibody (DPATB-H82199)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Goat Anti-Tobramycin Polyclonal Antibody
Target	Tobramycin
Immunogen	Tobramycin conjugated to BSA.
Isotype	IgG
Source/Host	Goat
Species Reactivity	N/A
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	ELISA
Format	Liquid
Size	200 µl
Buffer	Whole serum.
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C or -80°C (add glycerol to a final volume of 50% for extra stability). Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	Tobramycin is an aminoglycoside antibiotic used to treat various types of bacterial infections, particularly Gram-negative infections.
Keywords	Tobramycin; Anti – Tobramycin; 1-epitobramycin; 3'-deoxykanamycinb; d-6-tyrideoxy-alpha-d-ribohexopyranosyl-(1-6))-2-deoxy; distobram gernebcin lpha-d-ribohexopyranosyl-(1-4)]-2-deoxy-d-streptomine; nebcin; nebicina; nebramycin6; nebramycinvi