



# Anti-Trenbolone polyclonal antibody (DPATB-H82203)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	Sheep Anti-Trenbolone Polyclonal Antibody
<b>Target</b>	Trenbolone
<b>Immunogen</b>	Chemical / Small Molecule: 17 alpha Trenbolone-BTG.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Sheep
<b>Species Reactivity</b>	N/A
<b>Purification</b>	IgG fraction
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	cELISA
<b>Format</b>	Liquid
<b>Size</b>	50 µl
<b>Buffer</b>	20mM PBS, 150mM Sodium chloride, pH 7.2
<b>Preservative</b>	0.09% Sodium Azide
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

## BACKGROUND

<b>Introduction</b>	Trenbolone is a steroid used by veterinarians on livestock to increase muscle growth and appetite. To increase its effective half-life, trenbolone is not used in an unrefined form, but is rather administered as ester derivatives such as trenbolone acetate
<b>Keywords</b>	Trenbolone; 17β-Hydroxyestra-4,9,11-trien-3-one