



Anti-Gentamycin polyclonal antibody (DPATB-H82985)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Gentamycin Polyclonal Antibody
Target	Gentamycin
Immunogen	Chemical / Small Molecule: Gentamycin.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	RIA
Format	Liquid
Size	125 µl
Buffer	Whole serum.
Preservative	None
Storage	Store at +4°C short term. Aliquot and store at -20°C or -80°C (add glycerol to a final volume of 50% for extra stability). Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	Gentamycin is an aminoglycoside antibiotic used to treat infections caused by many different types of bacteria. Gentamycin is usually administered by intravenous infusion or intramuscular injection.
Keywords	Gentamicin;