



# Anti-Gentamycin polyclonal antibody (DPATB-H82986)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	Goat Anti-Gentamycin Polyclonal Antibody
<b>Target</b>	Gentamycin
<b>Immunogen</b>	Chemical / Small Molecule: Gentamicin.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Goat
<b>Species Reactivity</b>	N/A
<b>Purification</b>	Whole antiserum
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	RIA
<b>Format</b>	Liquid
<b>Size</b>	125 µl
<b>Buffer</b>	Whole serum.
<b>Preservative</b>	0.1% Sodium Azide
<b>Storage</b>	Add glycerol to a final volume of 50%, aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

## BACKGROUND

<b>Introduction</b>	Gentamycin is an aminoglycoside antibiotic used to treat infections caused by many different types of bacteria. Gentamycin is usually administered by intravenous infusion or intramuscular injection.
<b>Keywords</b>	Gentamicin;