



Anti-Taxol polyclonal antibody (DPATB-H83226)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Taxol Polyclonal Antibody
Target	Taxol
Immunogen	Taxol conjugation tg (Thyroglobulin bovine).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Bovine
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	ELISA
Format	Liquid
Size	250 µl
Buffer	Whole serum
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction Taxol (paclitaxel) is particularly effective as an anticancer agent because tubulin is not its only target. It also binds to a protein called Bcl2, which normally blocks the process of apoptosis, or cell death. Taxol therefore plays a dual role in destroy

Keywords Paclitaxel,