



Magic™ Anti-Raspberry bushy dwarf virus Coat Protein Polyclonal antibody (DPATB-H83588)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Chicken Anti-Raspberry Bushy Dwarf Virus Coat Protein Polyclonal Antibody
Target	Raspberry bushy dwarf virus Coat Protein
Immunogen	Native Raspberry Bushy Dwarf Virus (RBDV) particles purified from Chenopodium quinoa.
Isotype	IgY
Source/Host	Chicken
Species Reactivity	Raspberry bushy dwarf virus
Purification	IgY fraction
Conjugate	Unconjugated
Applications	WB, ELISA
Format	Liquid
Size	125 µg
Buffer	PBS, 1mM EDTA
Preservative	0.05% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	Raspberry bushy dwarf virus, which is spread by pollen and infects many red and black raspberry cultivars. It is also found in 'Boysenberry', 'Loganberry' and 'Marionberry' blackberries. The fruit from infected plants often will not make IQF grade and is
Keywords	RBDV;