



RHA™ anti-Streptomycin monoclonal antibody, clone STR (HMABPY045)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Cross reactivity: Dihydrostreptomycin 0.5ppb
Immunogen	Streptomycin [KLH]
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	STR
Purification	Purified with Absolute Mag™ Protein A Agarose (WHM-L086)
Conjugate	Unconjugated
Applications	ELISA, LFIA
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.