



RHA™ anti-Sudan monoclonal antibody, clone 112 (HMABPY081)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Cross-Reactivity: Para red 20ppb; Su-I 4ppb; Su-II 25ppb; Su-III 50ppb; Su-IV 50ppb
Immunogen	Sudan [KLH]
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	112
Purification	Purified with Absolute Mag™ Protein A Agarose (WHM-L086)
Conjugate	Unconjugated
Applications	ELISA, LFIA
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.